

THERA^{live}

PASSIVE AND ACTIVE EXERCISE FOR UPPER AND LOWER EXTREMITIES

You want to do at home what therapists recommend to improve your state of health?
You want to do daily exercise at home?

The new THERA-live gives you the opportunity to exercise whenever and wherever you want - even at home. Enjoy and manage your daily live easier by the effects of daily training.



FITS TO YOUR HOME
SMALL AND LIGHTWEIGHT



TOGETHER FOR A LIFE IN MOTION

» THERA-LIVE – SMALL, LIGHTWEIGHT AND EFFECTIVE.«

Who should use the THERA-live exerciser?

Integrating the new THERA-live into the everyday life is indicated at:

- Reduced mobility caused by multiple-sclerosis, para-/tetraplegia, stroke, Parkinson's disease
- Muscular diseases
- Orthopedic and rheumatic diseases
- Severe heart/blood circulation troubles
- Rehabilitation after hip/knee replacement

How does the new THERA-live help you?

Daily exercise with THERA-live Reduces spasticity. Therefore it eases things in daily life.

- Enhances strength and power to give you more stability and security within your mobility range.
- Activates your muscular power through our intelligent motor control.
- Reduces the consequences of a lack of movement such as cold legs, stiff joints, edemas
- Activates bladder and bowel system
- Affects your fitness and motivation.

Here are some testimonials:

MS patient:

»My knees and hips show more range of movement and my cold legs disappear.«

Paraplegic, L3:

»I can keep my muscular size much longer.«

MS patient:

»My spasticity is reduced for hours. So it makes me more independent.«

Parents of a head trauma patient:

»We are convinced that without the training his legs would be much more stiff.«

Person with Parkinson's disease:

»Walking is now much easier and movement is more smooth.«

Technical specifications:

- Control unit with display
- Easy START/STOP control
- Shows Watt, km, kcal, time, rpm, resistance level, activity of left/right leg
- Displays the activity of the left/right leg separately
- Spasticity control
- Integrated electronic flywheel
- Get-in-support
- 2-step adjustment pedal length
- Large foot cups
- Control unit holder
- Length: 54 cm, 21.3"
- Width: 46 cm, 18.1"
- Height: 46 cm, 18.1"



Enjoy your daily life. Supports your activities in daily life THERA-live integrates you into your family.



Easy to use. Ergonomically handles and buttons even for people with reduced motoric control.



Also available: Active/passive exercise for the upper tone and the arms. Please ask for our program of individual hand-grips and fixings.

Further options:
- Handle
- Calf-support
- Foot-fixation
- Tip up protection

How can I find out, if the THERA-live suits to my very individual needs?

The effects of THERA-live exercise may be very individual. Therefore we always give you the chance for free testing. Please feel free and give us a call. We will be happy to give you more information.



Unit 203, Milton Keynes Business Center, Foxhunter Drive
Linford Wood, Milton Keynes, MK14 6GD.
Tel: 01908 698895, Fax: 01908 698894,
Web: www.medicotech.co.uk, Email: info@medicotech.co.uk

